**Medical Engineering  
Technical Department**

**Letter to  
Med Manager Country  
Technical Manager Country**

Name            Geist  
Dept.            TD PS 1, (Med/LINA)  
Phone            +49 (9131) 84 3317  
Fax               +49 (9131) 84 7363  
E-Mail            PRODUCT TRACKING / LINA  
Date:             November 24, 2000

**Safety Modification SP017/00/S  
UROSKOP D1/D2/D3  
Safety Advisory Letter**

Reason for advisory:

Accidental System movements during operation of the UROSKOP D are reported (for example, the table tilts up, or the table top moves by itself).

Measures

Please distribute **immediately, by registered letter with reply coupon**, the letter „**Customer Safety Information – Activating the EMERGENCY STOP buttons in a hazardous situation**“ (see Enclosure) to the customers concerned.

Subsequently, please mark the field „MODIFICATION PERFORMED“ on the modification reply card and duly complete and return the latter to us **without delay**.

Equipment concerned:

Product:	Part No.	Type	Serial No.
UROSKOP D1	97 85 023	G5354	01002 – 01060
UROSKOP D2/D3	97 85 015	G5354	01001 – 01178

Material:

No material is required for this modification.

Very truly yours,

(Med/LINA)

Mailing Address:  
Siemens AG

TDPS1

P.O. Box 3260

Enclosure

### Customer Safety Information:

All Technical Service Departments with

UROSOP D installations (D1, D2, D3)

Name	Wilko Kuphal
Abteilung	SPU
Telefon	+49 91 31 84-37 34
Fax	+49 91 31 84-54 00
E-Mail	wilko.kuphal@med.siemens.de
Unser Zeichen	ku-fi
Datum	20. November 2000

#### re: Operating Notice for UROSOP D Systems (Activating the NOT –STOP Buttons in Hazardous Situations)

Dear Colleague:

As a result of customer complaints received here, we are asking you to draw the attention of customers operating UROSOP D systems in your sales region specifically to the fact that accidental movement of their UROSOP D system may occur during operation. You will find additional details in the sample letter, enclosed. You are responsible for translating this customer information into the local language and for informing the customers.

The reported possibility of error concerns a possible, sporadic error that may pose a danger to both the patient and the system.


The decision has been made to perform a software upgrade for all UROSOP D systems, starting in mid-December 2000. The corresponding Technical Service information is in preparation.

The enclosed customer letter is intended to notify the operators of UROSOP systems in advance and to ask them to be particularly attentive during UROSOP D operation.

The absolute rule is: **Whenever accidental system movement occurs (for example, table tilt or longitudinal table travel), system movement must immediately be interrupted by pressing an EMERGENCY STOP button, and Service must be notified.** (See Operating Instructions.)

Thank you for your understanding and for your active cooperation in informing our customers.

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### **Activating the EMERGENCY STOP buttons in a hazardous situation**

**Dear UROSKOP D User:**

Please observe the following important safety notice when operating your UROSKOP D system: In case of accidental system movement during operation of your UROSKOP D (for example, the system is moved from the horizontal into the upright position, or the tabletop moves by itself), immediately press one of the three red EMERGENCY STOP buttons on the system or the EMERGENCY SHUTDOWN switch in the room.

All system movements will be switched off immediately. If you have had to stop system movement in this way, please have the system inspected by our Customer Service personnel before resuming operation. (See also the Safety chapter in the Operating Instructions.)

We have now created an improved version of the upgrade software, and will begin installing this in mid-December 2000 at no charge on all UROSKOP D systems. Our Service personnel will be in contact with you to arrange a suitable date for implementation on your system.

Meanwhile, please instruct your operating staff of the need for particular attention when operating UROSKOP D.

Thank you very much for your understanding.

Siemens Aktiengesellschaft

Wilko Kuphal  
(SPU, Productmanager Urodiagnostic)

Ursula Sieberg  
(SPU, Marketingmanagement Urology)

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